

A Urinalysis System So Advanced It Lets You Do Less and See More

Atellica 1500 Automated Urinalysis System



Integrating the Power of Digital Imaging for Accurate, Reliable Urinalysis Testing

Automated urine chemistry and sediment testing in one streamlined unit


The Atellica™ 1500 Automated Urinalysis System* combines the trusted CLINITEK Novus® Automated Urine Chemistry Analyzer and the new Atellica™ UAS 800 Urine Sediment Analyzer in one completely automated unit. It's a new, streamlined system that lets you keep up and stay ahead—no matter how busy your urinalysis workload becomes.

The CLINITEK Novus analyzer combines the urine strip technology you have grown to trust with a cutting-edge digital camera to take your urine chemistry testing to new levels of clinical relevance.

The new Atellica UAS 800 analyzer is a completely automated urine sediment analyzer powered by advanced high-resolution digital imaging that closely replicates high-power field-like microscopy—designed to minimize the need for manual microscopy.

Helps to reduce everything but quality

The Atellica 1500 Automated Urinalysis System lets you manage more samples with less staff in a shorter time, while never compromising on high-quality results. Siemens Healthineers unique combination of urine chemistry and sediment technologies also streamlines the detection of UTIs and early kidney disease for better clinical outcomes.



Do less and see more
with a truly digital
automated urinalysis
system that sets the new
standard for accuracy
and efficiency.

**The Atellica 1500 Automated Urinalysis System:
Just what you'd expect from
the pioneers in modern urinalysis.**

Atellica 1500 Automated Urinalysis System*			
CLINITEK Novus Analyzer with CLINITEK Novus Cassettes		Atellica UAS 800 Analyzer*	
Menu	CLINITEK Novus 10 Cassette <ul style="list-style-type: none"> • Bilirubin • Clarity • Color • Glucose • Ketone • Leukocyte esterase • Nitrite • Occult blood • pH • Protein • Specific gravity • Urobilinogen 	CLINITEK Novus PRO 12 Cassette* <ul style="list-style-type: none"> • Albumin • Albumin-to-creatinine ratio • Bilirubin • Clarity • Color • Creatinine • Glucose • Ketone • Leukocyte esterase • Nitrite • Occult blood • pH • Protein • Protein-to-creatinine ratio • Specific gravity • Urobilinogen 	Auto-reported Categories 11 Major Categories: <ul style="list-style-type: none"> • Bacteria • Crystals • Hyaline cast • Mucus • Nonsquamous EC • Pathological cast • RBC • Sperm • Squamous EC • WBC • Yeast Sub-categories: <ul style="list-style-type: none"> • Bacteria cocci • Bacteria rods • WBC clumps
	Measurement Technology	Digital camera that measures the color changes in the entire RGB color spectrum and applies color correction for improved accuracy	

*Not available for sale in the U.S. Product availability varies by country.

Atellica Portfolio of Laboratory Products

Engineered by Siemens Healthineers to deliver control and simplicity so you can drive better outcomes.

Tighter control of your lab, simplified workflow, and more time to focus on driving better business and clinical outcomes—that's the promise of our Atellica™ portfolio of laboratory products.

Control.
Simplicity.
Better Outcomes.

From the point of care to the central lab, Siemens Healthineers has performance-driven urinalysis solutions to meet your needs.

The same reagent strip technology standardizes testing at every location. Providing early detection, monitoring, and management of diseases, Siemens Healthineers offers the widest range of urinalysis products on the market.

More than
75 Years
of Urinalysis
Leadership and Innovation

Atellica, CLINITEK Novus, and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc., or its affiliates. All other trademarks and brands are the property of their respective owners.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone: +49 9131 84 0
siemens.com/healthineers

Local Contact Information

Siemens Healthcare Diagnostics Inc.
Point of Care Diagnostics
2 Edgewater Drive
Norwood, MA 02062-4637
USA
Phone: +781 269-3000