

Atellica UAS 60 Analyzer

Powerful Imaging, Now in a Smaller Size

Transform your urine testing workflow
with accelerated sediment analysis

siemens-healthineers.com



SIEMENS
Healthineers

Empower Your Lab with a Faster Workflow

Atellica UAS 60 Analyzer

Accelerate your urine testing workflow with the Atellica® UAS 60 Analyzer. This new system features the same innovative technology used in the Atellica UAS 800 Analyzer in a smaller footprint. The Atellica UAS 60 Analyzer can improve the quality and consistency of your urine sediment results to help revolutionize your lab's capabilities while saving time, cost, and space.

Trusted technology. Consistent results.

Powered by advanced high-resolution digital imaging, you can automatically analyze your patient samples with less effort and more accuracy. The Atellica UAS 60 Analyzer produces results that closely replicate the quality of manual microscopic analysis. It uses high-quality imaging to generate up to 15 full-field images, effectively reducing the time it takes to analyze and process samples compared to the traditional manual method.

Maximize efficiency with minimal training

- ✓ Automated urine sediment analysis in a small footprint
- ✓ Fast time to result compared to manual microscopy: up to 60 samples per hour
- ✓ No liquid reagents: requires only a 3-minute cuvette column change
- ✓ Advanced algorithm to perform particle counting and result reporting for highly accurate urine sediment analysis
- ✓ Full-field-of-view digital images for simple review: no manual microscopy required



Make informed diagnostic decisions with industry-leading digital imaging


Urine sediment testing as easy as 1, 2, 3



Initialize the testing process with a few simple steps that require minimal training.



Process samples automatically after pipetting the sample onto the cuvette.



Perform automated particle counting and results reporting with optional manual viewing and identification mode.

Atellica UAS 60 Analyzer* Menu

<p>Auto-reported categories</p>	<p>Major categories</p> <ul style="list-style-type: none"> • Bacteria • Crystals • Hyaline cast • Mucus • Nonsquamous EC • Pathological cast • Squamous EC • WBC • Yeast 	<p>Subcategories</p> <ul style="list-style-type: none"> • Bacteria cocci • Bacteria rods • WBC clumps <p>Crystal subcategories</p> <ul style="list-style-type: none"> • CaOxm • CaOxd • URI • TRI
<p>Measurement technology</p>	<ul style="list-style-type: none"> • Digital microscopy imaging • Dual-focusing mechanism to generate clear images • Real-time image capture function 	

See More with Less Effort

With high-quality digital imaging, deliver consistent results with confidence.

*Not available for sale in the U.S. Product availability varies by country.

Legacy of Innovation

Siemens Healthineers offers one of the widest ranges of urinalysis products on the market. From the point of care to the central lab, Siemens Healthineers has performance-driven urinalysis solutions to meet your needs. We are continuously advancing industry-leading technology to produce confidence-inspiring results for early detection, monitoring, and management of diseases.

Family of Urine Analyzers



Atellica 1500 Automated Urinalysis System†

Atellica UAS 60 Analyzer

**More than 80 Years of
Urinalysis Leadership
and Innovation**



We are continuously advancing
industry-leading technology for
more confidence in results.

†The Atellica 1500 Automated Urinalysis System is a combination of the Atellica UAS 800 Analyzer and CLINITEK Novus® Analyzer.

At Siemens Healthineers, we pioneer breakthroughs in healthcare. For everyone. Everywhere. By constantly bringing breakthrough innovations to market, we enable healthcare professionals to deliver high-quality care, leading to the best possible outcome for patients.

Our portfolio, spanning from in-vitro and in-vivo diagnostics to image-guided therapy and innovative cancer care, is crucial for clinical decision-making and treatment pathways. With our strengths in patient twinning, precision therapy, as well as digital, data, and artificial intelligence (AI), we are well positioned to take on the biggest challenges in healthcare. We will continue to build on these strengths to help fight the world's most threatening diseases, improving the quality of outcomes, and enabling access to care.

We are a team of 66,000 highly dedicated employees across more than 70 countries passionately pushing the boundaries of what's possible in healthcare to help improve people's lives around the world.

Atellica, CLINITEK, and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc., or its affiliates. All other trademarks and brands are the property of their respective owners.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone: +49 9131 84-0
siemens-healthineers.com

Published by

Siemens Healthcare Diagnostics Inc.
Point of Care Diagnostics
2 Edgewater Drive
Norwood, MA 02062-4637
USA
Phone: +1 781-269-3000